





Donnenfeld

Femto Instruments



Watch it!



Donnenfeld Femto Spatula

This new spatula is designed to open femtosecond created cataract incisions. Its 1mm long, flattened, "spear-shaped" tip aids in identifying and entering the interface and opening the wound with a sweeping motion. This instrument easily opens incisions which cannot be opened with traditional femto spatulas.



Donnenfeld Femto Splitter

This splitter is designed specially to compliment modern cataract femtosecond ablation techniques. The paracentral curve in the shaft exactly matches the shape of a laser created central nuclear bowl providing optimal surface contact for efficiently splitting the lens out to the periphery and down to the posterior plate. The Donnenfeld Femto Splitter improves nuclear chopping with minimal manipulation.

Designed by
Eric D. Donnenfeld, MD
of Rockville Centre, NY

 **katena**
DESIGNED FOR SIGHT®

973-989-1600 • 800-225-1195 • www.katena.com